I SEMESTER PRACTICAL MANUAL FOR FIRST B.Sc. CHEMISTRY (w. e. f. 2020 – 2021)

PRACTICAL – 1(ANALYSIS OF SALT MIXTURE)



PREPARED BY S. ANIL DEV M. Sc. LECTURER IN CHEMISTRY

DEPARTMENT OF CHEMISTRY D.N.R.COLLEGE (AUTONOMOUS) BHIMAVARAM

SYLLABUS

LABORATORY COURSE -I 30hrs (2 H / W)

Practical-I Analysis of SALT MIXTURE (At the end of Semester-I)

Qualitative inorganic analysis (Minimum of Six mixtures should be analyzed)50 MCourse outcomes:

At the end of the course, the student will be able to;

- Understand the basic concepts of qualitative analysis of inorganic mixture
- Use glassware, equipment and chemicals and follow experimental procedures in the laboratory
- Apply the concepts of common ion effect, solubility product and concepts related to qualitative analysis

ANALYSIS OF SALT MIXTURE

50 M

Analysis of mixture salt containing two anions and two cations (From two different groups) from the following:

Anions

Carbonate, Sulphate, Chloride, Bromide, Acetate, Nitrate, Borate, Phosphate.

Cations

Lead, Copper, Iron, Aluminium, Zinc, Nickel, Manganese, Calcium, Strontium, Barium, Potassium and Ammonium.

SCHEME OF VALUATION

Time: 3 Hours	Maximum Marks: 50 M
Record	Marks:5M
Viva-Voce	Marks: 5M
Practical	Marks: 40M
For Preliminary Examination	Marks: 8M
1. Colour	1 Mark
2. Structure	1 Mark
3. Solubility	2 Marks
4. Action of Heat	2 Marks
5. Flame colour test	2 Marks
For anion	Marks: 16M
6. Identification test for anion	4 Marks
7. Confirmation tests for anion without soda	6 Marks
(For phosphate and sulphate 6 marks shall be added to the	marks at item 8 without this test)

8.	Procedure for the preparation of soda extract	2 Marks
9.	Confirmatory test for anion with soda extract	4 Marks

(For carbonate and borate 4 marks shall be added to the marks at item 6 without this test)

'or cation	Marks:14M
10. Identification of the Group	4 Marks
11. Colour of the ppt. in the identified Group	2 Marks
12. Confirmatory test for cation in the individual Group	8 Marks
	ACCARD IN THE READ AND CARD ON

Report

2 Marks

S. N. A. HOD OF CHEMISTRY Head of the Dept. of Chemistry. D.N.R. COLLEGE BHIMAVARAM-534 202.

QUALITATIVE INORGANIC SALT MIXTURE ANALYSIS

Experiment	Observation	Inference
	White	Salt mixtures of Cu, Fe, Ni and Mn are absent
	Blue or bluish green	May be a copper salt mixture
Colour	Brown or yellowish brown	May be a ferric salt mixture
	Pale green	May be a salt mixture of ferrous or nickel
	Green	May be CuCO ₃ or CuCl ₂
	Pale pink	May be a manganese salt mixture
Structure	Crystalline	May be a hydrated salt mixture
	White amorphous (appears as face powder)	May be ZnCO ₃ , BaCO ₃ , SrCO ₃ CaCO ₃
	Green amorphous	May be CuCO ₃
	White deliquescent salt mixture (wet with moisture in air)	May be ZnCl ₂ , Zn(NO ₃) ₂ , Ba(NO ₃) ₂ or CaCl ₂
	Pink deliquescent salt mixture	May be MnCl ₂
	Water drops formed on cooler parts of the test tube	Hydrated salt mixture
Action of heat	White coat formed on sides of the test tube (sublimation took place) and ammonia smell evolved.	May be an ammonium halide
	ernel salar annon	distant individent data

	White sublimate formed Without ammonia smell.	May be an aluminium salt mixture
	Colourless vapours evolved with vinegar smell.	May be acetate
	Light brown vapours evolved with pungent smell.	May be nitrate
(M	Yellow in cold, brown in hot	May be a lead salt mixture
	Blue in cold, white in hot.	May be copper sulphate
	Light yellow in cold, black in hot.	May be ferric chloride
an.	White in cold, orange yellow in hot.	May be a zinc salt mixture
	Green in cold, yellow in hot.	May be a nickel salt mixture
Flame colouration test On a watch glass a little of	Livid blue (dark blue)	Lead
the salt mixture was mixed with few drops of conc. HCl	Bluish green	Copper
The paste was exposed to	Green flashes	Zinc
colourless flame by a platinum wire or glass rod.	Apple green(light green)	Barium
observed which indicated	Crimson red	Strontium
presence of the cation	Brick red	Calcium
	Violet	Potassium

IDENTIFICATION TESTS FOR ANIONS

1. Test with dil. HCl: To a little of the salt mixture dil. HCl was added	Colourless and odourless gas (CO ₂) with brisk Effervescence evolved. On passing the gas lime water became milky white (CaCO ₃).	Carbonate (CO ₃ 2) present
	<u>On heating the test tube</u> smell of vinegar (CH ₃ COOH) observed.	Acetate (CH ₃ COO ⁻) present
2. Oxalic acid rubbing test for acetate: On a watch glass a little of the salt mixture, oxalic acid and a few drops of water were taken and rubbed with a glass rod or thumb.	Vinegar smell observed.	A cetate present
3. Test with conc. H₂SO4: In a dry test tube, to a little of the salt mixture a few drops of conc. H ₂ SO ₄ was added. (without heating)	White fumes evolved or Colourless pungent gas (HCl) evolved, fumed in moist air (or when blown by air from mouth)	Chloride (Cl [®]) present
	Reddish brown vapours evolved	Bromide (Br) present
With heating	vinegar smell observed	Acetate present
with heating	Light reddish brown vapours evolved.	Nitrate (NO ₃) present
4. Test for Borate: On a watch glass, to a little of the salt mixture the same quantity of Calcium fluoride (CaF ₂) powder was mixed and made a paste by adding a few drops of conc. H ₂ SO ₄	Green edged flame observed. (the volatile BF ₃ formed in the test imparted green colour to flame)	Borate (BO ₃ ³⁻ or B ₂ O ₇ ^{2–}) present

and exposed the paste to colourless flame with the help of a glass rod.		
5. Test for Phosphate: A little of the salt mixture was dissolved in a minimum of conc. HNO ₃ by heating, and excess of ammonium molybdate added.	A canary yellow precipitate, (NH4)3PO4.12MoO3, formed.	Phosphate (PO ₄ ^{3–}) present
6. Test for Sulphate: To the salt mixture solution barium chloride solution was added.	A white precipitate (BaSO ₄) formed, insoluble in conc. HCl or conc. HNO ₃	Sulphate (SO4 ^{2—}) present

CONFIRMATORY TESTS FOR ANIONS

1. BaCl ₂ test for carbonate: To the salt mixture solution (prepared in water) BaCl ₂ solution was added.	White precipitate (BaCO ₃) formed, soluble in dil. HCl	Carbonate confirmed
2. Ethyl or Amyl acetate test for acetate: In a dry test tubea little of the salt mixture, a few drops of conc. H ₂ SO ₄ and ethyl alcohol or amyl alcohol were taken and warmed. The contents in the test tube were poured in a beaker full of water.	A fruity smell of ethyl acetate (CH ₃ COOC ₂ H ₅) observed. or Ripen banana smell of amyl acetate (CH ₃ COOC ₅ H ₁₁) observed.	Acetate confirmed
3. MnO ₂ test for halides: In a dry test tube, a little of the saltmixture and manganese dioxide were wet with a few drops of	Greenish yellow chlorine gas (Cl ₂) evolved. Reddish brown bromine gas (Br ₂) evolved.	Chloride confirmed Bromide confirmed
conc. H_2SO_4 and then heated.		

<u>4. Chromyl chloride test</u> <u>for chloride:</u> In a dry test tube a little of the salt mixture and potassium dichromate were taken and wet with a few drops of conc. H ₂ SO ₄ and then it was heated.	Dark red chromyl chloride vapours (CrO_2Cl_2) formed and passed downwards from test tube mouth. Nearby the mouth another test tube containing lead acetate solution was brought. The vapours were dissolved in the solution and yellow precipitate (PbCrO ₄) formed.	Chloride confirmed
5. Copper turnings test for nitrate: In a dry test tube a little of the salt mixture, few copper turnings and few drops of conc. H_2SO_4 were taken and heated.	Reddish brown vapours evolved. or Formation of light brown gas (NO ₂) intensified. Mass in the test tube turned green (due to formation of anhydrous copper nitrate).	Nitrate (NO ₃ -) present
6. Ethyl borate test for borate: In a dry test tube a little of the salt mixture, few drops of conc. H_2SO_4 and ethyl alcohol were taken and warmed. The vapours evolved were exposed to flame and continued warming of the test tube.	The vapours were continued to burn with green edge flame (the volatile ethyl borate formed imparted green colour to the flame).	Borate confirmed

Soda Extract tests for further confirmation

Preparation of sodium carbonate extract:

The salt mixture was mixed with twice or thrice the amount of sodium carbonate. Sufficient water was added, boiled and filtered. The filtrate is known as sodium carbonate extract or **soda extract**.

<u>1.Test for sulphate:</u> A small portion of the soda extract was taken. On adding dil. HCl formation of bubbles started. Addition of the acid continued until the bubbles stopped (neutralisation). BaCl ₂ solution added.	White precipitate (BaSO ₄) formed insoluble in conc. HCl or conc. HNO ₃ .	Sulphate (SO4 ²⁻) confirmed
2. Test for acetate : The soda extract was neutralized with dil. HCl and neutral ferric chloride solution was added. (Add NaHCO ₃ to FeCl ₃ solution until bubbles stopped to prepare neutral FeCl ₃ solution)	Red precipitate, Fe(CH₃COO) ₃ , formed.	Acetate confirmed
3. Confirmatory tests for halides: A small portion of the soda extract was neutralized with dil. HNO ₃ . Then AgNO ₃ solution added.	White precipitate (AgCl) formed, completely soluble in ammonia solution Pale yellow precipitate (AgBr) formed, sparingly soluble in ammonia solution	Chloride confirmed Bromide confirmed
4. Brown ring test for nitrate: A small portion of the soda extract was neutralized with dil. H ₂ SO ₄ . Then freshly prepared FeSO ₄ solution was added in the same quantity. To this Conc. H ₂ SO ₄ was added drop by drop from sides of the test tube by keeping it in inclined position.	A brown ring, [Fe(NO)]SO4, formed at the junction of the two layers.	Nitrate confirmed

TESTS FOR CATION

Preparation of Original Solution

The solvent at top of the list in observation column for the solubility test was taken. The given salt mixture was added to it little by little and shaken well to get a clear solution. It is 0.S.

<u>Tests for ammonium cation (NH₄⁺)</u>

Experiment	Observation	Inference
1. NaOH test: To a little of the salt mixture NaOH was solution added and warmed.	Colourless gas with ammonia smell (NH ₃ gas) observed. A glass rod dipped in HCl solution when exposed to this gas gave dense white fumes (NH ₄ Cl).	Ammonium ion present
2.Nessler's reagent test : To the salt mixture solution excess of NaOH solution and few drops of Nessler's reagent (alkaline solution of K ₂ HgI ₄) were added.	Reddish-brown (chocolate colour) precipitate or solution, NH ₂ .HgO.HgI.H ₂ O (Iodide of Millon's base), formed.	Ammonium ion confirmed

GROUP SEPARATION TABLE

To t	he original :	solution dil.HCl i	is added and	centrifuged		
Residue:	Centrifuga	ate: H_2S gas is pa	assed throug	h the centrif	ugate and centrifuged	
PbCl ₂ (white ppt)	Residue: Black – CuS Brown - BiS Yellow - CdS	Centrifugate: ' added and cer	Гhe centrifu ıtrifuged	gate is boiled	l off till H_2S is eliminated. After cooling excess of NH_4Cl and NH_4OH are	
		- Residue: - White gelatanious- F Al(OH) ₃ Dirty green- W Fe(OH) ₂ Brown- F Fe(OH) ₃ N E N	Centrifugate: H Residue: Cen of N White -	ete: H ₂ S gas i Centrifuga of NH ₄ Cl, N	H ₂ S gas is passed through the centrifugate and centrifuged. entrifugate: Centrifugate is boiled off till H ₂ S is eliminated. After cooling excess f NH ₄ Cl, NH ₄ OH and (NH ₄) ₂ CO ₃ are added and centrifuged.	
			Flesh - MnS Black – NiS,CoS	Residue: White- BaCO ₃ SrCO ₃ CaCO ₃	Centrifugate: After boiling, this is used in the identification of VI group cations.	

NAME OF THE FACULTY: S. ANIL DEV M. Sc.

LECTURER IN CHEMISTRY D.N.R.COLLEGE (A), BHIMAVARAM

I GROUP ANALYSIS

To the ppt. obtained in I group distilled water was added, boiled and filtered.

Residue	Filtrate				
is not					
formed	It was divided into two parts.				
	Part I: To this part potassium chromate solution was added. Yellow ppt. (PbCrO ₄) formed. Lead ion (Pb²⁺) present.				
	Part II: (Golden spangles test: To this part potassium iodide solution was added. Yellow ppt. (PbI) formed. To the separated ppt. distilled water was added, boiled, and on cooling under tap water, particles like golden yellow spangles (PbI) were appeared. Lead ion (Pb²⁺) present .				

II GROUP ANALYSIS

(Pb²⁺ and Cu²⁺ of Group IIA cations present)

The ppt. obtained in II group was washed with water, 33% nitric acid (1 part conc. acid + 2 parts water) was added, heated for a minute, dil. sulphuric acid was added and filtered.

Residue		Filtrate		
Ammonium acetate solution was added, heated and filtered.		Ammonia solution was added until ammonia smells, warmed and filtered.		
Residue	Filtrate	Residue	Filtrate	
Mercuric mercury Hg ²⁺	Acidified with dilute acetic acid and potassium chromate solution was added. Yellow ppt. (PbCrO ₄) formed. Lead ion (Pb ²⁺) present.	Bismuth ion Bi ³⁺	To this blue colour solution acetic acid and potassium ferrocyanide solution were added. Reddish brown ppt. (Chocolate colour ppt.) Cu ₂ [Fe(CN) ₆] formed. Cupric ion (Cu²⁺) present.	

III GROUP ANALYSIS

To the ppt. obtained in III group NaOH solution was added in excess and filtered.

Residue	Filtrate	
Dil. HCl was added to dissolve it and potassium ferrocyanide solution was	Dil. HCl was added, and ammonia solution was added until ammonia smells.	
ppt. or solution, Fe ₄ [Fe(CN) ₆] _{3,} formed.	White gelatinous ppt., Al(OH) ₃ , formed.	
Iron ion (Fe ²⁺ or Fe ³⁺) present.	Aluminium ion (Al ³⁺) present.	

IV GROUP ANALYSIS

To the ppt. obtained in IV group dil. HCl was added, warmed and filtered.

Residue	Filtrate		
It was dissolved in minimum of aqua regia (3 parts Conc. HCl + 1 part Conc. HNO ₃) and was diluted by water	It was boiled to expel H ₂ S gas. NaOH solution was added in excess, boiled and filtered. Residue Filtrate		
Test for Nickel: To the above solution ammonia solution was added until ammonia smells. Few drops of Di Methyl Glyoxime (DMG) reagent was added. Rose coloured ppt., [Ni (DMG) ₂] complex, formed.	Conc.HNO ₃ was added to dissolve the residue, lead dioxide was added, boiled, cooled and diluted with water. Purple pink colour, Mn(OH) ₂ , formed.	H ₂ S was passed. White ppt. (ZnS) formed.	
Nickel ion, Ni ²⁺ present.	Manganese ion (Mn ²⁺) present.	Zinc ion (Zn ²⁺) present.	

V GROUP ANALYSIS

The ppt. obtained in V group was dissolved in dil. acetic acid and was divided into three parts.

Part I	Part II	Part III
Potassium chromate solution was added. Yellow ppt. (BaCrO ₄) formed.	Ammonia solution was added until ammonia smells, and ammonium sulphate solution was added and shaken well. White ppt. (SrSO ₄) formed.	Ammonia solution was added until ammonia smells, and ammonium oxalate solution was added. White ppt. (CaC ₂ O ₄) formed.
Barium ion (Ba ²⁺) present	Strontium ion (Sr ²⁺) present.	Calcium ion (Ca²+) present.

VI GROUP ANALYSIS

Before going to the General Group Table ammonium ion was tested. So, tests were done for the remaining cation in this VI group potassium cation (K⁺) only.

Experiment	Observation	Inference
1.To the O.S. acetic acid and sodium cobaltinitrite solution were added.	Yellow ppt. K ₃ [Co(NO ₂) ₆] formed.	Potassium ion (K+) present.
2.To the O.S. tartaric acid solution was added and sides of the test tube were scratched with a glass rod.	White ppt., (COOH)(CHOH)₂(COOK), formed	Potassium ion (K+) present.

Report: The given salt mixture contains the following ions.

Anions:

Cations: